

# Venosan Legline Size Chart

## Measurement Guidelines

In order to work effectively graduated compression hosiery must be the correct size and apply the right amount of pressure to the leg. It is essential that accurate measurements are taken.

## Circumference


- All measurements should be taken on the bare leg, in centimetres, whilst standing (if possible).
- Ensure that the tape measure is snug but not too tight.
- Measurements should be taken when the leg is rested and not swollen

## Leg Length

This product has Maximum leg lengths for each style of hose. Please refer to measurements guidelines and chart below for instructions.

## Measurement Points

1. The length of the leg from the base of the heel to just below the knee **ID**
2. Ankle - Ankle at the narrowest point above the ankle bone **cB**
3. Calf - Widest point at the calf **cC**
4. Below Knee - Narrowest point below the knee **cD**
5. Mid Thigh - The circumference of the leg in the middle of the thigh **cF**



## Circumference Chart

Circumference (cm)	X-Small	Small	Medium	Large	X-large
cB	20-21cm 7.9-8.3"	22-23cm 8.7-9.1"	24-25cm 9.4-9.8"	26-27cm 10.2-10.6"	28-29cm 11-11.4"
cC	31-33cm 12.2-13"	34-36cm 13.4-14.2"	37-39cm 14.6-15.4"	40-42cm 15.7-16.5"	43-45cm 16.9-17.7"
cD	29-31cm 11.4-12.2"	32-34cm 12.6-13.4"	35-37cm 13.8-14.6"	38-40cm 15-15.7"	41-43cm 16.1-16.9"
cF	42-44cm 16.5-17.3"	45-47cm 17.7-18.5"	48-50cm 18.9-19.7"	51-53cm 20.1-20.9"	54-60cm 21.3-23.6"